

For Research Use Only

GLUD1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67026-6-PBS



Basic Information

Catalog Number: 67026-6-PBS	GenBank Accession Number: BC040132	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 2746	CloneNo.: 1E1C9
Source: Mouse	UNIPROT ID: P00367	
Isotype: IgG2b	Full Name: glutamate dehydrogenase 1	
Immunogen Catalog Number: AG6179	Calculated MW: 61 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50458-3, GLUT1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67026-6-PBS. Detection antibody: 67026-5-PBS. Standard: Ag6179. Range: 0.781-100 ng/mL.