

For Research Use Only

# HAI-1 Monoclonal antibody, PBS Only



Catalog Number: 67031-1-PBS

## Basic Information

Catalog Number:

67031-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG26733

GenBank Accession Number:

BC004140

GeneID (NCBI):

6692

UNIPROT ID:

O43278

Full Name:

serine peptidase inhibitor, Kunitz type

1

Calculated MW:

58 kDa

Observed MW:

58 kDa

Purification Method:

Protein A purification

CloneNo.:

2H2D8

## Applications

Tested Applications:

WB, Indirect ELISA, IF

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



human placenta tissue were subjected to SDS PAGE followed by western blot with 67031-1-Ig (HAI-1 antibody) at dilution of 1:900 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67031-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (-20°C Methanol) fixed MCF-7 cells using HAI-1 antibody (67031-1-Ig, Clone: 2H2D8) at dilution of 1:800 and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L). This data was developed using the same antibody clone with 67031-1-PBS in a different storage buffer formulation.