

For Research Use Only

Integrin alpha-4 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67040-2-PBS



Basic Information

Catalog Number: 67040-2-PBS	GenBank Accession Number: BC156712	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 3676	CloneNo.: 1G4G11
Source: Mouse	UNIPROT ID: P13612	
Isotype: IgG2a	Full Name: integrin, alpha 4 (antigen CD49D, alpha 4 subunit of VLA-4 receptor)	
Immunogen Catalog Number: AG28024	Calculated MW: 1032 aa, 115 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50496-1, Integrin alpha-4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67040-2-PBS. Detection antibody: 67040-3-PBS. Standard: Ag28024. Range: 1.563-50 ng/mL.