

For Research Use Only

# Integrin alpha-4 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67040-3-PBS



## Basic Information

Catalog Number:

67040-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG28024

GenBank Accession Number:

BC156712

GeneID (NCBI):

3676

UNIPROT ID:

P13612

Full Name:

integrin, alpha 4 (antigen CD49D,  
alpha 4 subunit of VLA-4 receptor)

Calculated MW:

1032 aa, 115 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1E12B10

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50496-1, Integrin alpha-4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67040-2-PBS. Detection antibody: 67040-3-PBS. Standard: Ag28024. Range: 1.563-50 ng/mL.