

For Research Use Only

SPATA7 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67046-1-PBS



Basic Information

Catalog Number: 67046-1-PBS	GenBank Accession Number: BC008656	Purification Method: Protein G purification
Size: 1mg/ml	GeneID (NCBI): 55812	CloneNo.: 1F6F1
Source: Mouse	UNIPROT ID: Q9POW8	
Isotype: IgG1	Full Name: spermatogenesis associated 7	
Immunogen Catalog Number: AG21234	Calculated MW: 599 aa, 68 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51085-3, SPATA7 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67046-2-PBS. Detection antibody: 67046-1-PBS. Standard:Ag21234. Range: 0.391-100 ng/mL.