

For Research Use Only

IL4R Monoclonal antibody, PBS Only (Detector)



Catalog Number: 67051-2-PBS

Basic Information

| | | |
|---|---|---|
| Catalog Number: 67051-2-PBS | GenBank Accession Number: NM_000418 | Purification Method: Protein G purification |
| Size: 1 mg/ml | GeneID (NCBI): 3566 | CloneNo.: 1E9B5 |
| Source: Mouse | UNIPROT ID: P24394 | |
| Isotype: IgG1 | Full Name: interleukin 4 receptor | |
| Immunogen Catalog Number: AG28411 | Calculated MW: 90 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50242-1, IL4R Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67051-3-PBS. Detection antibody: 67051-2-PBS. Standard:Ag28411. Range: 1.563-100 ng/mL.



Cytometric bead array standard curve of MP50242-4, IL4R Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67051-4-PBS. Detection antibody: 67051-2-PBS. Standard:Ag28411. Range: 0.391-25 ng/mL.