

For Research Use Only

# IL4R Monoclonal antibody, PBS Only (Capture/Detector)



Catalog Number: 67051-3-PBS

## Basic Information

Catalog Number:

67051-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgM

Immunogen Catalog Number:

AG28411

GenBank Accession Number:

NM\_000418

GeneID (NCBI):

3566

UNIPROT ID:

P24394

Full Name:

interleukin 4 receptor

Calculated MW:

90 kDa

Purification Method:

Euglobulin precipitation

CloneNo.:

3F12E8

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


## Selected Validation Data



Cytometric bead array standard curve of MP50242-1, IL4R Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67051-3-PBS. Detection antibody: 67051-2-PBS. Standard:Ag28411. Range: 1.563-100 ng/mL



Cytometric bead array standard curve of MP50242-2, IL4R Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67051-4-PBS. Detection antibody: 67051-3-PBS. Standard:Ag28411. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP50242-3, IL4R Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67051-3-PBS. Detection antibody: 67051-5-PBS. Standard:Ag28411. Range: 0.195-100 ng/mL