

For Research Use Only

# Tissue Factor Monoclonal antibody, PBS Only



Catalog Number: 67056-1-PBS

## Basic Information

Catalog Number:

67056-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG27563

GenBank Accession Number:

BC011029

GeneID (NCBI):

2152

UNIPROT ID:

P13726

Full Name:

coagulation factor III (thromboplastin,  
tissue factor)

Calculated MW:

295 aa, 33 kDa

Observed MW:

45-50 kDa

Purification Method:

Protein A purification

CloneNo.:

2E4A1

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



A431 cells were subjected to SDS PAGE followed by western blot with 67056-1-Ig (Tissue Factor antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67056-1-PBS in a different storage buffer formulation.