

For Research Use Only

F3 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67056-3-PBS



Basic Information

Catalog Number:

67056-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgM

Immunogen Catalog Number:

AG27563

GenBank Accession Number:

BC011029

GeneID (NCBI):

2152

UNIPROT ID:

P13726

Full Name:

coagulation factor III (thromboplastin,

tissue factor)

Calculated MW:

295 aa, 33 kDa

Purification Method:

Euglobulin precipitation

CloneNo.:

3G11A7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50510-1, F3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67056-2-PBS. Detection antibody: 67056-3-PBS. Standard: Ag27563. Range: 0.195-25 ng/mL