

For Research Use Only

F3 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67056-4-PBS



Basic Information

Catalog Number: 67056-4-PBS	GenBank Accession Number: BC011029	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 2152	CloneNo.: 4A9C7
Source: Mouse	UNIPROT ID: P13726	
Isotype: IgG1	Full Name: coagulation factor III (thromboplastin, tissue factor)	
Immunogen Catalog Number: AG27563	Calculated MW: 295 aa, 33 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50510-2, F3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67056-4-PBS. Detection antibody: 67056-5-PBS. Standard: Ag27563. Range: 1.563-100 ng/mL.