

For Research Use Only

REDD1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67059-3-PBS



Basic Information

Catalog Number:

67059-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG0965

GenBank Accession Number:

BC007714

GeneID (NCBI):

54541

UNIPROT ID:

Q9NX09

Full Name:

DNA-damage-inducible transcript 4

Calculated MW:

25 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2C4A6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50649-1, REDD1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67059-2-PBS. Detection antibody: 67059-3-PBS. Standard:Ag0965. Range: 0.098-50 ng/mL



Cytometric bead array standard curve of MP50649-2, REDD1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67059-4-PBS. Detection antibody: 67059-3-PBS. Standard:Ag0965. Range: 0.098-50 ng/mL