

For Research Use Only

REDD1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67059-6-PBS



Basic Information

Catalog Number:

67059-6-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG0965

GenBank Accession Number:

BC007714

GeneID (NCBI):

54541

UNIPROT ID:

Q9NX09

Full Name:

DNA-damage-inducible transcript 4

Calculated MW:

25 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1H5G6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50649-3, REDD1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67059-5-PBS. Detection antibody: 67059-6-PBS. Standard: Ag0965. Range: 0.098-50 ng/mL