

For Research Use Only

Myogenin Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67082-3-PBS



Basic Information

Catalog Number: 67082-3-PBS	GenBank Accession Number: BC053899	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 4656	CloneNo.: 4E10B5
Source: Mouse	UNIPROT ID: P15173	
Isotype: IgG1	Full Name: myogenin (myogenic factor 4)	
Immunogen Catalog Number: AG6414	Calculated MW: 25 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51032-1, Myogenin Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67082-2-PBS. Detection antibody: 67082-3-PBS. Standard: Ag6414. Range: 0.781-50 ng/mL.



Cytometric bead array standard curve of MP51032-2, Myogenin Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67082-1-PBS. Detection antibody: 67082-3-PBS. Standard: Ag6414. Range: 1.563-25 ng/mL.