

For Research Use Only

RALB Monoclonal antibody, PBS Only (Capture)



Catalog Number: 67094-2-PBS

Basic Information

Catalog Number: 67094-2-PBS	GenBank Accession Number: BC018163	Purification Method: Protein G purification
Size: 1mg/ml	GeneID (NCBI): 5899	CloneNo.: 1G7F2
Source: Mouse	UNIPROT ID: P11234	Affinity: $K_D = 1.10 \times 10^{-9}$
Isotype: IgG1	Full Name: v-ral simian leukemia viral oncogene homolog B (ras related; GTP binding protein)	$K_{off} = 1.89 \times 10^{-5}$ $K_{on} = 1.72 \times 10^4$
Immunogen Catalog Number: AG28600	Calculated MW: 206 aa, 23 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50129-1, RALB Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67094-2-PBS. Detection antibody: 67094-1-PBS. Standard:Ag28600. Range: 6.25-100 ng/mL