

For Research Use Only

DUSP12 Monoclonal antibody, PBS Only



Catalog Number: 67101-1-PBS

Basic Information

Catalog Number: 67101-1-PBS	GenBank Accession Number: BC006286	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 11266	CloneNo.: 1G8G10
Source: Mouse	UNIPROT ID: Q9UNI6	
Isotype: IgG1	Full Name: dual specificity phosphatase 12	
Immunogen Catalog Number: AG8241	Calculated MW: 370 aa, 43 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HEK-293 cells were subjected to SDS PAGE followed by western blot with 67101-1-Ig (DUSP12 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67101-1-PBS in a different storage buffer formulation.