

For Research Use Only

CYP3A4 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67110-4-PBS



Basic Information

Catalog Number: 67110-4-PBS	GenBank Accession Number: BC101631	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 1576	CloneNo.: 2E6E5
Source: Mouse	UNIPROT ID: P08684	
Isotype: IgG3	Full Name: cytochrome P450, family 3, subfamily A, polypeptide 4	
Immunogen Catalog Number: AG13074	Calculated MW: 503 aa, 57 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50718-2, CYP3A4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67110-2-PBS. Detection antibody: 67110-4-PBS. Standard: Ag13074. Range: 0.195-100 ng/mL



Cytometric bead array standard curve of MP50718-3, CYP3A4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67110-5-PBS. Detection antibody: 67110-4-PBS. Standard: Ag13074. Range: 0.781-100 ng/mL