

For Research Use Only

CYP3A4 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67110-5-PBS



Basic Information

Catalog Number:

67110-5-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG13074

GenBank Accession Number:

BC101631

GeneID (NCBI):

1576

UNIPROT ID:

P08684

Full Name:

cytochrome P450, family 3, subfamily A, polypeptide 4

Calculated MW:

503 aa, 57 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1G7B1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50718-3, CYP3A4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67110-5-PBS. Detection antibody: 67110-4-PBS. Standard: Ag13074. Range: 0.781-100 ng/mL.