

For Research Use Only

PKC Delta Monoclonal antibody, PBS Only



Catalog Number: 67113-1-PBS

Basic Information

Catalog Number: 67113-1-PBS	GenBank Accession Number: BC043350	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 5580	CloneNo.: 2B3E3
Source: Mouse	UNIPROT ID: Q05655	
Isotype: IgG2b	Full Name: protein kinase C, delta	
Immunogen Catalog Number: AG5631	Calculated MW: 78 kDa	
	Observed MW: 78-80 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human, mouse, rat

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 67113-1-Ig (PKC delta antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67113-1-PBS in a different storage buffer formulation.