

For Research Use Only

SPRY2 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67123-1-PBS



Basic Information

Catalog Number: 67123-1-PBS	GenBank Accession Number: BC015745	Purification Method: Protein G purification
Size: 1mg/ml	GeneID (NCBI): 10253	CloneNo.: 1E12C1
Source: Mouse	UNIPROT ID: O43597	
Isotype: IgG1	Full Name: sprouty homolog 2 (Drosophila)	
Immunogen Catalog Number: AG28599	Calculated MW: 35 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human, mouse, rat

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50527-1, SPRY2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67123-1-PBS. Detection antibody: 67123-2-PBS. Standard: Ag28599. Range: 0.391-50 ng/mL