

For Research Use Only

# DDR2 Monoclonal antibody, PBS Only



Catalog Number: 67126-1-PBS

## Basic Information

<b>Catalog Number:</b> 67126-1-PBS	<b>GenBank Accession Number:</b> NM_001014796	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 4921	<b>CloneNo.:</b> 1E9A10
<b>Source:</b> Mouse	<b>UNIPROT ID:</b> Q16832	
<b>Isotype:</b> IgG1	<b>Full Name:</b> discoidin domain receptor tyrosine kinase 2	
<b>Immunogen Catalog Number:</b> AG28340	<b>Calculated MW:</b> 97 kDa	
	<b>Observed MW:</b> 100-110 kDa	

## Applications

**Tested Applications:**  
WB, ELISA

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Jurkat cells were subjected to SDS PAGE followed by western blot with 67126-1-Ig (DDR2 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67126-1-PBS in a different storage buffer formulation.