

For Research Use Only

ABCA7 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67128-3-PBS



Basic Information

Catalog Number:

67128-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG27499

GenBank Accession Number:

NM_019112

GeneID (NCBI):

10347

UNIPROT ID:

Q8IZY2

Full Name:

ATP-binding cassette, sub-family A
(ABC1), member 7

Calculated MW:

234 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2E12A2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51162-1, ABCA7 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67128-2-PBS. Detection antibody: 67128-3-PBS. Standard:Ag27499. Range: 3.125-200 ng/mL.