

For Research Use Only

STAR Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 67130-3-PBS



Basic Information

Catalog Number:

67130-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG28698

GenBank Accession Number:

BC010550

GeneID (NCBI):

6770

UNIPROT ID:

P49675

Full Name:

steroidogenic acute regulatory
protein

Calculated MW:

285 aa, 32 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1D5B1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP50772-1, STAR Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67130-2-PBS. Detection antibody: 67130-3-PBS. Standard:Ag28698. Range: 0.195-200 ng/mL



Cytometric bead array standard curve of MP50772-2, STAR Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67130-3-PBS. Detection antibody: 67130-4-PBS. Standard:Ag28698. Range: 0.195-200 ng/mL



Cytometric bead array standard curve of MP50772-3, STAR Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67130-3-PBS. Detection antibody: 67130-5-PBS. Standard:Ag28698. Range: 6.25-200 ng/mL