

For Research Use Only

Ninein Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67132-3-PBS



Basic Information

Catalog Number: 67132-3-PBS	GenBank Accession Number: BC034708	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 51199	CloneNo.: 1E5A5
Source: Mouse	UNIPROT ID: Q8N4C6	
Isotype: IgG1	Full Name: ninein (GSK3B interacting protein)	
Immunogen Catalog Number: AG14652	Calculated MW: 2090 aa, 243 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50497-1, Ninein Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67132-2-PBS. Detection antibody: 67132-3-PBS. Standard: Ag14652. Range: 0.098-100 ng/mL.