

For Research Use Only

TRIM21 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 67136-3-PBS

Basic Information

Catalog Number: 67136-3-PBS	GenBank Accession Number: BC010861	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 6737	CloneNo.: 1D11D10
Source: Mouse	UNIPROT ID: P19474	
Isotype: IgG1	Full Name: tripartite motif-containing 21	
Immunogen Catalog Number: AG28377	Calculated MW: 475 aa, 54 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50371-1, TRIM21 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67136-2-PBS. Detection antibody: 67136-3-PBS. Standard:Ag28377. Range: 0.098-50 ng/mL



Cytometric bead array standard curve of MP50371-2, TRIM21 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67136-1-PBS. Detection antibody: 67136-3-PBS. Standard:Ag28377. Range: 0.781-100 ng/mL