

For Research Use Only

COQ6 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67162-2-PBS



Basic Information

Catalog Number: 67162-2-PBS	GenBank Accession Number: BC014181	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 51004	CloneNo.: 2D9E2
Source: Mouse	UNIPROT ID: Q9Y2Z9	
Isotype: IgG1	Full Name: coenzyme Q6 homolog, monooxygenase (S. cerevisiae)	
Immunogen Catalog Number: AG28991	Calculated MW: 468 aa, 51 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50616-1, COQ6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67162-1-PBS. Detection antibody: 67162-2-PBS. Standard: Ag28991. Range: 0.098-100 ng/mL.