

For Research Use Only

# COQ6 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67162-4-PBS



## Basic Information

<b>Catalog Number:</b> 67162-4-PBS	<b>GenBank Accession Number:</b> BC014181	<b>Purification Method:</b> Protein A Magarose purification
<b>Size:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 51004	<b>CloneNo.:</b> 3C1H10
<b>Source:</b> Mouse	<b>UNIPROT ID:</b> Q9Y2Z9	
<b>Isotype:</b> IgG2a	<b>Full Name:</b> coenzyme Q6 homolog, monooxygenase (S. cerevisiae)	
<b>Immunogen Catalog Number:</b> AG28991	<b>Calculated MW:</b> 468 aa, 51 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage


**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:  
T: 4006900926      E: Proteintech-CN@ptglab.com      W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50616-2, COQ6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67162-3-PBS. Detection antibody: 67162-4-PBS. Standard:Ag28991. Range: 1.563-100 ng/mL.