

For Research Use Only

TRIB2 Monoclonal antibody

Catalog Number: 67171-1-Ig **1 Publications**



Basic Information

Catalog Number:

67171-1-Ig

Size:

1700 µg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG7501

GenBank Accession Number:

BC002637

GeneID (NCBI):

28951

UNIPROT ID:

Q92519

Full Name:

tribbles homolog 2 (Drosophila)

Calculated MW:

39 kDa

Observed MW:

35-40 kDa

Purification Method:

Protein G purification

CloneNo.:

1A10A6

Recommended Dilutions:

WB 1:1000-1:6000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

CoIP

Species Specificity:

Human, Mouse

Cited Species:

human

Positive Controls:

WB : A549 cells, U-937 cells, RAW 264.7 cells

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Sam A Jamieson	36108635	Structure	CoIP

Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



A549 cells were subjected to SDS PAGE followed by western blot with 67171-1-Ig (TRIB2 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.