

For Research Use Only

F11R Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67186-1-PBS



Basic Information

Catalog Number: 67186-1-PBS	GenBank Accession Number: BC001533	Purification Method: Protein G purification
Size: 1mg/ml	GeneID (NCBI): 50848	CloneNo.: 1F12F10
Source: Mouse	UNIPROT ID: Q9Y624	
Isotype: IgG1	Full Name: F11 receptor	
Immunogen Catalog Number: AG25656	Calculated MW: 33 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50597-1, F11R Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67186-1-PBS. Detection antibody: 67186-2-PBS. Standard: Ag25656. Range: 0.195-50 ng/mL