

For Research Use Only

SHC1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67194-4-PBS



Basic Information

Catalog Number:

67194-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG25173

GenBank Accession Number:

NM_183001

GeneID (NCBI):

6464

UNIPROT ID:

P29353

Full Name:

SHC (Src homology 2 domain
containing) transforming protein 1

Calculated MW:

63 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1H4G2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50763-2, SHC1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67194-3-PBS. Detection antibody: 67194-4-PBS. Standard: Ag25173. Range: 0.391-100 ng/mL.



Cytometric bead array standard curve of MP50763-3, SHC1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67194-2-PBS. Detection antibody: 67194-4-PBS. Standard: Ag25173. Range: 0.391-50 ng/mL.