

For Research Use Only

SHC1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67194-5-PBS



Basic Information

Catalog Number: 67194-5-PBS	GenBank Accession Number: NM_183001	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 6464	CloneNo.: 2C11A4
Source: Mouse	UNIPROT ID: P29353	
Isotype: IgG2a	Full Name: SHC (Src homology 2 domain containing) transforming protein 1	
Immunogen Catalog Number: AG25173	Calculated MW: 63 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50763-4, SHC1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67194-5-PBS. Detection antibody: 67194-3-PBS. Standard: Ag25173. Range: 0.781-100 ng/mL.