

For Research Use Only

EIF4G1 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 67199-2-PBS

Basic Information

Catalog Number: 67199-2-PBS	GenBank Accession Number: BC007788	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 1981	CloneNo.: 1C3G8
Source: Mouse	UNIPROT ID: Q04637	
Isotype: IgG2b	Full Name: eukaryotic translation initiation factor 4 gamma, 1	
Immunogen Catalog Number: AG8357	Calculated MW: 1600 aa, 176 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50357-1, EIF4G1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67199-2-PBS. Detection antibody: 67199-3-PBS. Standard: Ag8357. Range: 6.25-100 ng/mL