

For Research Use Only

TMEM111 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67205-3-PBS



Basic Information

Catalog Number:

67205-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG18365

GenBank Accession Number:

BC022807

GeneID (NCBI):

55831

UNIPROT ID:

Q9POL2

Full Name:

transmembrane protein 111

Calculated MW:

261 aa, 30 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2B11A4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51179-1, TMEM111 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67205-2-PBS. Detection antibody: 67205-3-PBS. Standard:Ag18365. Range: 0.098-100 ng/mL



Cytometric bead array standard curve of MP51179-2, TMEM111 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67205-4-PBS. Detection antibody: 67205-3-PBS. Standard:Ag18365. Range: 0.098-100 ng/mL