

For Research Use Only

# FGFR1OP Monoclonal antibody, PBS Only (Capture)

Catalog Number:67207-2-PBS



## Basic Information

Catalog Number: 67207-2-PBS	GenBank Accession Number: BC011902	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 11116	CloneNo.: 1C12E7
Source: Mouse	UNIPROT ID: O95684	
Isotype: IgG1	Full Name: FGFR1 oncogene partner	
Immunogen Catalog Number: AG28671	Calculated MW: 43 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50541-1, FGFR1OP Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67207-2-PBS. Detection antibody: 67207-3-PBS. Standard: Ag28671. Range: 0.781-100 ng/mL