

For Research Use Only

FGFR1OP Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67207-3-PBS



Basic Information

Catalog Number: 67207-3-PBS	GenBank Accession Number: BC011902	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 11116	CloneNo.: 1E10A9
Source: Mouse	UNIPROT ID: O95684	
Isotype: IgG2b	Full Name: FGFR1 oncogene partner	
Immunogen Catalog Number: AG28671	Calculated MW: 43 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50541-1, FGFR1OP Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67207-2-PBS. Detection antibody: 67207-3-PBS. Standard: Ag28671. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP50541-2, FGFR1OP Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67207-4-PBS. Detection antibody: 67207-3-PBS. Standard: Ag28671. Range: 0.781-100 ng/mL.