

For Research Use Only

FGFR1OP Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67207-5-PBS



Basic Information

| | | |
|---|--|--|
| Catalog Number: 67207-5-PBS | GenBank Accession Number: BC011902 | Purification Method: Protein A Magarose purification |
| Size: 1 mg/ml | GeneID (NCBI): 11116 | CloneNo.: 1G8A2 |
| Source: Mouse | UNIPROT ID: O95684 | |
| Isotype: IgG2b | Full Name: FGFR1 oncogene partner | |
| Immunogen Catalog Number: AG28671 | Calculated MW: 43 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50541-3, FGFR1OP Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67207-4-PBS. Detection antibody: 67207-5-PBS. Standard: Ag28671. Range: 1.563-100 ng/mL