

For Research Use Only

DTX2 Monoclonal antibody, PBS Only



Catalog Number: 67209-1-PBS

Basic Information

Catalog Number: 67209-1-PBS	GenBank Accession Number: BC064535	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 113878	CloneNo.: 2D4F8
Source: Mouse	UNIPROT ID: Q86UW9	
Isotype: IgG2a	Full Name: deltex homolog 2 (Drosophila)	
Immunogen Catalog Number: AG12022	Calculated MW: 253aa, 28 kDa; 575aa, 64 kDa	
	Observed MW: 68 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67209-1-Ig (DTX2 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67209-1-PBS in a different storage buffer formulation.