

For Research Use Only

# TMEM230 Monoclonal antibody

Catalog Number: 67247-1-Ig



## Basic Information

Catalog Number:

67247-1-Ig

Size:

1500 µg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG15080

GenBank Accession Number:

BC011990

GeneID (NCBI):

29058

UNIPROT ID:

Q96A57

Full Name:

chromosome 20 open reading frame

30

Calculated MW:

183 aa, 20 kDa

Observed MW:

16 kDa

Purification Method:

Protein A purification

CloneNo.:

4D3A7

Recommended Dilutions:

WB 1:5000-1:50000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human, Mouse, Rat

Positive Controls:

WB : HeLa cells, COLO 320 cells, Caco-2 cells, HT-29 cells, A549 cells, MCF-7 cells, HSC-T6 cells, NIH/3T3 cells, PC-12 cells, HEK-293 cells

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67247-1-Ig (TMEM230 antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours.