

For Research Use Only

B4GALT3 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67254-3-PBS



Basic Information

Catalog Number:

67254-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgM

Immunogen Catalog Number:

AG28038

GenBank Accession Number:

BC009985

GeneID (NCBI):

8703

UNIPROT ID:

O60512

Full Name:

UDP-Gal:betaGlcNAc beta 1,4-galactosyltransferase, polypeptide 3

Calculated MW:

44 kDa

Purification Method:

Euglobulin precipitation

CloneNo.:

2B6B1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51163-1, B4GALT3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67254-2-PBS. Detection antibody: 67254-3-PBS. Standard: Ag28038. Range: 0.391-100 ng/mL



Cytometric bead array standard curve of MP51163-2, B4GALT3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67254-4-PBS. Detection antibody: 67254-3-PBS. Standard: Ag28038. Range: 1.563-50 ng/mL