

For Research Use Only

SIRT5 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67257-2-PBS



Basic Information

Catalog Number:

67257-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG7743

GenBank Accession Number:

BC000126

GeneID (NCBI):

23408

UNIPROT ID:

Q9NXA8

Full Name:

sirtuin (silent mating type information regulation 2 homolog) 5 (S. cerevisiae)

Calculated MW:

34 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2D3G3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50563-1, SIRT5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67257-2-PBS. Detection antibody: 67257-3-PBS. Standard: Ag7743. Range: 0.098-100 ng/mL.