

For Research Use Only

PSG1 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 67269-2-PBS

Basic Information

Catalog Number: 67269-2-PBS	GenBank Accession Number: BC022338	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 5669	CloneNo.: 1C3B2
Source: Mouse	UNIPROT ID: P11464	
Isotype: IgG1	Full Name: pregnancy specific beta-1- glycoprotein 1	
Immunogen Catalog Number: AG28604	Calculated MW: 426 aa, 48 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50302-1, PSG1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67269-2-PBS. Detection antibody: 67269-3-PBS. Standard: Ag28604. Range: 0.781-100 ng/mL.