

For Research Use Only

PSG1 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 67269-3-PBS

Basic Information

Catalog Number:

67269-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG28604

GenBank Accession Number:

BC022338

GeneID (NCBI):

5669

UNIPROT ID:

P11464

Full Name:

pregnancy specific beta-1-
glycoprotein 1

Calculated MW:

426 aa, 48 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2E11B3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50302-1, PSG1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67269-2-PBS. Detection antibody: 67269-3-PBS. Standard: Ag28604. Range: 0.781-100 ng/mL.