

For Research Use Only

# CD5 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 67278-3-PBS

## Basic Information

<b>Catalog Number:</b> 67278-3-PBS	<b>GenBank Accession Number:</b> BC027901	<b>Purification Method:</b> Protein G purification
<b>Size:</b> 1mg/ml	<b>GeneID (NCBI):</b> 921	<b>CloneNo.:</b> 1A10A5
<b>Source:</b> Mouse	<b>ENSEMBL Gene ID:</b> ENSG00000110448	
<b>Isotype:</b> IgG1	<b>UNIPROT ID:</b> P06127	
<b>Immunogen Catalog Number:</b> AG11097	<b>Full Name:</b> CD5 molecule	
	<b>Calculated MW:</b> 495 aa, 55 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.

**Storage Buffer:**  
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50028-1, CD5 Monoclonal Matched Antibody Pair - BSA and Azide free. Capture antibody: 67278-2-PBS. Detection antibody: 67278-3-PBS. Standard: Ag11097. Range: 12.5-100 ng/mL.