

For Research Use Only

ALG13 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67283-2-PBS



Basic Information

Catalog Number:

67283-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG14955

GenBank Accession Number:

BC005336

GeneID (NCBI):

79868

UNIPROT ID:

Q9NP73

Full Name:

asparagine-linked glycosylation 13
homolog

Calculated MW:

165 aa, 18 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1E3E3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50876-1, ALG13 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67283-1-PBS. Detection antibody: 67283-2-PBS. Standard: Ag14955. Range: 3.125-200 ng/mL.