

For Research Use Only

# Rat Insulin1 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 67284-3-PBS

## Basic Information

Catalog Number:

67284-3-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG28809

GenBank Accession Number:

NM\_019129.3

GeneID (NCBI):

24505

UNIPROT ID:

P01322

Full Name:

insulin 1

Calculated MW:

12 kDa

Purification Method:

Protein A purification

CloneNo.:

1D11C2

## Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

Rat

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Sandwich ELISA standard curve of MP50081-1, Rat Insulin1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67284-2-PBS. Detection antibody: 67284-3-PBS. Standard: Ag28809. Range: 0.313-20 ng/mL