

For Research Use Only

WFS1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67291-2-PBS



Basic Information

Catalog Number: 67291-2-PBS	GenBank Accession Number: BC030130	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 7466	CloneNo.: 3D11G5
Source: Mouse	UNIPROT ID: O76024	
Isotype: IgG3	Full Name: Wolfram syndrome 1 (wolframin)	
Immunogen Catalog Number: AG25735	Calculated MW: 890 aa, 100 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50883-1, WFS1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60620-1-PBS. Detection antibody: 67291-2-PBS. Standard:Ag2114. Range: 6.25-100 ng/mL