

For Research Use Only

luciferase Monoclonal antibody, PBS Only



Catalog Number: 67293-1-PBS

Basic Information

Catalog Number:

67293-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG13000

GenBank Accession Number:

AY015993

GeneID (NCBI):

Full Name:

Calculated MW:

20 kDa

Purification Method:

Protein A purification

CloneNo.:

2D2E1

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Gaussia princeps

Background Information

This antibody is designed to detect the luciferase of luciferase (Gaussia princeps) tagged protein.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Luciferase tagged protein were subjected to SDS PAGE followed by western blot with 67293-1-Ig (Luciferase antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67293-1-PBS in a different storage buffer formulation.