

For Research Use Only

KLF6 Monoclonal antibody, PBS Only



Catalog Number: 67297-1-PBS

Basic Information

Catalog Number:

67297-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG6859

GenBank Accession Number:

BC000311

GeneID (NCBI):

1316

UNIPROT ID:

Q99612

Full Name:

Kruppel-like factor 6

Calculated MW:

32 kDa

Observed MW:

32-42 kDa

Purification Method:

Protein A purification

CloneNo.:

2F8B3

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human, mouse

Background Information

KLF6 (kruppel-like factor 6) is a zinc finger transcription factor and tumor suppressor with biological activities and transcriptional targets in growing range. It is highly expressed in placenta followed by spleen, thymus, prostate, testis, small intestine and colon. One key target gene of KLF6 is endoglin, which is a homodimeric cell membrane glycoprotein and TGF- β auxiliary receptor, owns a pro-angiogenic role in endothelial cells and is also involved in malignant progression. Also KLF6 induces apoptosis in prostate cancer cells through upregulation of ATF3. KLF6 was normally detected as a 46 kDa protein, which is larger than translation product (32 kDa), because of post-translational modification including both phosphorylation and glycosylation. This antibody raises against the full length of KLF6 gene of human origin, and can recognize both of KLF6, including KLF6 (~32 kDa), p-KLF6 (~35-46 kDa). Otherwise, there may be a non-specific band (50 kDa) shown in detection.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67297-1-Ig (KLF6 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67297-1-PBS in a different storage buffer formulation.