

For Research Use Only

MYH1 Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 67299-2-PBS



Basic Information

Catalog Number:

67299-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG17129

GenBank Accession Number:

BC114545

GeneID (NCBI):

4619

UNIPROT ID:

P12882

Full Name:

myosin, heavy chain 1, skeletal muscle, adult

Calculated MW:

1939 aa, 223 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2D12G4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50449-1, MYH1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67299-2-PBS. Detection antibody: 67299-3-PBS. Standard: Ag17129. Range: 1.563-100 ng/mL.



Cytometric bead array standard curve of MP50449-3, MYH1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67299-6-PBS. Detection antibody: 67299-2-PBS. Standard: Ag17129. Range: 0.391-25 ng/mL.