

For Research Use Only

# MYH1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67299-5-PBS



## Basic Information

|   |  |   |
|---|--|---|
| <b>Catalog Number:</b><br>67299-5-PBS       | <b>GenBank Accession Number:</b><br>BC114545                       | <b>Purification Method:</b><br>Protein G purification |
| <b>Size:</b><br>1 mg/ml                     | <b>GeneID (NCBI):</b><br>4619                                      | <b>CloneNo.:</b><br>1E5C4                             |
| <b>Source:</b><br>Mouse                     | <b>UNIPROT ID:</b><br>P12882                                       |   |
| <b>Isotype:</b><br>IgG1                     | <b>Full Name:</b><br>myosin, heavy chain 1, skeletal muscle, adult |   |
| <b>Immunogen Catalog Number:</b><br>AG17129 | <b>Calculated MW:</b><br>1939 aa, 223 kDa                          |   |

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50449-2, MYH1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67299-4-PBS. Detection antibody: 67299-5-PBS. Standard: Ag17129. Range: 3.125-100 ng/mL.