

For Research Use Only

TLR8 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67317-3-PBS



Basic Information

Catalog Number: 67317-3-PBS	GenBank Accession Number: BC101075	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 51311	CloneNo.: 2B9A5
Source: Mouse	UNIPROT ID: Q9NR97	
Isotype: IgG1	Full Name: toll-like receptor 8	
Immunogen Catalog Number: AG16566	Calculated MW: 1041 aa, 120 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50542-1, TLR8 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67317-2-PBS. Detection antibody: 67317-3-PBS. Standard: Ag16566. Range: 0.098-100 ng/mL.



Cytometric bead array standard curve of MP50542-2, TLR8 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67317-4-PBS. Detection antibody: 67317-3-PBS. Standard: Ag16566. Range: 0.098-100 ng/mL.